

DuraSonic .75 Gallon Digital Ultrasonic Cleaner | DS3L



Stock Photo - Options May Vary

FACTORY NEW

DuraSonic DS3L .75 Gallon Digital Ultrasonic Cleaner. This UltraSonic Cleaner is an Industrial Grade Model with a Sweep and Degas Function. Temperature Controlled Heat for best cleaning performance. Additional features include Operator Temperature Controls (knob), Digital Ultrasonic Timer, and a very Reliable Industrial-Grade Transducer. Do you have questions?

To view a wide variety of Autoclaves & Equipment, click the following link to browse from automatic autoclaves and distillers to manual autoclaves and ultrasonic cleaners. Venture Medical's Complete Catalog of Autoclaves.

FEATURES

- Industrial Grade Model
- Sweep and Degas Function
- Temperature Controlled Heat for best cleaning performance
- Operator Temperature Controls (knob)
- Digital Ultrasonic Timer

SPECIFICATIONS

- Tank Capacity: .75 gal / 2.8 ltr
- Inside Dimensions (LxWxH): 9"x5.5"x4"
- Outside Dimensions (L x W x H): 12"x 7"x 10.5"
- Sonic Power Watt Average: 120

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)