

# DuraSonic .75 Gallon Ultrasonic Cleaner w/Heat & Mech Timer | DS230HT



Stock Photo - Options May Vary

## FACTORY NEW

DuraSonic DS230HT .75 Gallon Ultrasonic Cleaner w/Heat and Mechanical Timer. The new benchmark in benchtop cleaning!

Durasonic is the Sound Choice and is Now Available with Ceramically Enhanced Transducers and Lifetime. DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a benchtop price.

The driving force is the power source: industrial stack transducers offer superior reliability, but benchtop cleaners use wafer-thin discs. Durasonic brings full-wave 45 kHz Ceramically Enhanced Industrial Transducers to the Benchtop line. All units feature 0-30 minute analog or 0-99 minute digital timers.

Additionally, thermostatically controlled heaters improve efficiency. When comparing benchtop cleaners, focus on features—there is no comparison.

## FEATURES

- Analog model with timer and adjustable thermostatically controlled heater
- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency 45 kHz
- Adjustable, thermostatically controlled heaters (ambient to 80°C)

## SPECIFICATIONS

- Capacity: .75 Gallons
- Dimensions: 9.5" x 5.25" x 4"
- All Stainless Steel Enclosures
- Lifetime Heater Warranty
- 2-year warranty on parts and labor

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)