

DuraSonic 1 Gallon Digital Ultrasonic Cleaner | DS360D



Stock Photo - Options May Vary

FACTORY NEW

DuraSonic DS360D 1 Gallon Digital Ultrasonic Cleaner. The new benchmark in benchtop cleaning! Durasonic is the Sound Choice and is Now Available with Ceramically Enhanced Transducers and Lifetime. DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a benchtop price.

The driving force is the power source. While industrial stack transducers offer superior reliability, benchtop cleaners use wafer-thin discs for ultrasonic energy. Durasonic brings full-wave 45 kHz Ceramically Enhanced Industrial Transducers to the benchtop line. All units feature 0-30 minute analog or 0-99 minute digital timers.

Additionally, thermostatically controlled heaters enhance efficiency. When comparing benchtop cleaners, focus on features—there is no comparison.

FEATURES

- Digital model with timer, heat, power control, and degas
- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency 45 kHz
- Adjustable, thermostatically controlled heaters (ambient to 80°C)

SPECIFICATIONS

- Capacity: 1 Gallon
- Dimensions: 9.5" x 5.25" x 6"
- All Stainless Steel Enclosures
- Timer: 0 - 99 Minutes
- Lifetime Heater Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)