

# DuraSonic 1 Gallon Ultrasonic Cleaner w/Heat & Mech Timer | DS360HT



Stock Photo - Options May Vary

## FACTORY NEW

DuraSonic DS360HT Ultrasonic Cleaner w/Heat & Mech Timer. The new benchmark in benchtop cleaning! DuraSonic is the Sound Choice and is Now Available with Ceramically Enhanced Transducers and Lifetime.

DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a benchtop price. The driving force is the power source: While industrial stack transducers are known for reliability, benchtop cleaners rely on wafer-thin discs. DuraSonic brings the superior performance of full-wave 45 kHz Ceramically Enhanced Industrial Transducers to the Benchtop line.

All units include 0-30 minute analog or 0-99 minute digital timers. Additionally, the thermostatically controlled heaters enhance efficiency. When comparing benchtop cleaners, focus on features—there is no comparison.

## FEATURES

- Analog model with timer and adjustable thermostatically controlled heater
- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency 45 kHz
- Adjustable, thermostatically controlled heaters (ambient to 80°C)

## SPECIFICATIONS

- Capacity: 1 Gallon
- Dimensions: 9.5" x 5.25" x 6"
- All Stainless Steel Enclosures
- Timer: 0 - 30 Minutes
- Lifetime Heater Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)