

DuraSonic 1.5 Gal Digital Ultrasonic Cleaner | DS6L



Stock Photo - Options May Vary

FACTORY NEW

DuraSonic DS6L 1.5 Gal Digital Ultrasonic Cleaner. The DS6L is an Industrial Grade Model Sweep and Degas Function with Temperature Controlled Heat for best cleaning performance. Features include an Operator Temperature Controls (knob) along with a Digital Ultrasonic Timer. DuraSonic Ultrasonic Cleaners have a very Reliable Industrial-Grade Transducer. Do you have questions?

To view a wide variety of Autoclaves & Equipment, click the following link to browse automatic autoclaves and distillers to manual autoclaves and ultrasonic cleaners. Venture Medical's Complete Catalog of Autoclaves.

FEATURES

- DuraSonic Digital Ultrasonic Cleaner - 1.5 Gallon (1.5 liters)
- Industrial Grade Model
- Sweep and Degas Function
- Temperature Controlled Heat for best cleaning performance
- Operator Temperature Controls (knob)

SPECIFICATIONS

- Inside Dimensions (L x W x H): 11.75" x 6" x 6"
- Outside Dimensions (L x W x H): 14" x 7" x 12.5"
- Sonic Power Watt Average: 180
- Includes Heater: Yes

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)