

DuraSonic 1.5 Gallon Digital Ultrasonic Cleaner | DS500D



Stock Photo - Options May Vary

FACTORY NEW

DuraSonic DS500D 1.5 Gallon Digital Ultrasonic Cleaner. The new benchmark in benchtop cleaning! Durasonic is the Sound Choice and is Now Available with Ceramically Enhanced Transducers and Lifetime.

DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a benchtop price. The driving force is the power source: While industrial stack transducers provide reliability, benchtop cleaners rely on wafer-thin discs for ultrasonic energy. Additionally, DuraSonic brings full-wave 45 kHz Ceramically Enhanced Industrial Transducers to the benchtop line.

All units include either 0-30 minute analog or 0-99 minute digital timers. The thermostatically controlled heaters enhance efficiency. When comparing benchtop cleaners, focus on features—there is no comparison.

FEATURES

- Digital model with timer, heat, power control, and degas
- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency 45 kHz
- Adjustable, thermostatically controlled heaters (ambient to 80°C)

SPECIFICATIONS

- Capacity: 1.5 Gallons
- Dimensions: 11.75" x 6" x 6"
- All Stainless Steel Enclosures
- Timer: 0 - 99 Minutes
- Lifetime Heater Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)