

# DuraSonic 2.5 Gallon Digital Ultrasonic Cleaner | DS1200D



Stock Photo - Options May Vary

## FACTORY NEW

DuraSonic DS1200D 2.5 Gallon Digital Ultrasonic Cleaner. The new benchmark in benchtop cleaning! DuraSonic stands as the sound choice, now available with Ceramically Enhanced Transducers and a lifetime warranty. DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a benchtop price.

The driving force is the power source: Industrial stack transducers are reliable, but benchtop cleaners use wafer-thin discs for energy. DuraSonic brings superior performance with full-wave 45 kHz Ceramically Enhanced Industrial Transducers to benchtop cleaners. Additionally, units come with either a 0-30 minute analog or 0-99 minute digital timer.

Heaters are thermostatically controlled for optimal efficiency.

## FEATURES

- Digital model with timer, heat, power control, and degas
- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency 45 kHz
- Adjustable, thermostatically controlled heaters (ambient to 80°C)

## SPECIFICATIONS

- Capacity: 2.5 Gallons
- Dimensions: 19.5" x 5.25" x 6"
- All Stainless Steel Enclosures
- Timer: 0 - 99 Minutes
- Lifetime Heater Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)