

DuraSonic 2.5 Gallon Ultrasonic Cleaner w/Heat & Mech Timer | DS1200HT



Stock Photo - Options May Vary

FACTORY NEW

DuraSonic DS1200HT 2.5 Gallon Ultrasonic Cleaner w/Heat & Mech Timer. The new benchmark in benchtop cleaning! DuraSonic is the Sound Choice and is Now Available with Ceramically Enhanced Transducers and Lifetime. DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a benchtop price.

The driving force is the power source: Industrial stack transducers excel in reliability and performance. Benchtop cleaners use wafer-thin discs. DuraSonic brings full-wave 45 kHz Ceramically Enhanced Industrial Transducers to the Benchtop line. Each unit includes a 0-30 minute analog or 0-99 minute digital timer.

The thermostatically controlled heaters enhance efficiency. When comparing benchtop cleaners, focus on features—there is no comparison.

FEATURES

- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency 45 kHz
- Adjustable, thermostatically controlled heaters (ambient to 80°C)
- Drains with a flow control valve

SPECIFICATIONS

- Capacity: 2.5 Gallons
- Dimensions: 19.5" x 5.25" x 6"
- All Stainless Steel Enclosures
- Lifetime Heater Warranty
- 2-year warranty on parts and labor

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)