

DuraSonic 3.25 Gallon Ultrasonic Cleaner w/Heat & Mech. Timer | DS1100HT



Stock Photo - Options May Vary

FACTORY NEW

DuraSonic DS1100HT 3.25 Gallon Ultrasonic Cleaner w/Heat and Mechanical Timer. The new benchmark in benchtop cleaning! DuraSonic is the sound choice, now available with Ceramically Enhanced Transducers and a lifetime warranty. DuraSonic's Rugged Industrial-Strength Performance gives you consistent quality at a benchtop price.

The driving force is the power source: Industrial stack transducers deliver superior reliability, while benchtop cleaners use wafer-thin discs. DuraSonic now brings full-wave 45 kHz Ceramically Enhanced Industrial Transducers to benchtop models. Each unit includes either a 0-30 minute analog or 0-99 minute digital timer.

The heaters are thermostatically controlled for efficiency. When comparing benchtop cleaners, compare features—no other unit compares to DuraSonic's performance.

FEATURES

- Analog model with timer and adjustable thermostatically controlled heater
- Ceramically Enhanced Transducers
- All stainless steel enclosures resist damage from chemical spills and will not crack
- Sweep Frequency 45 kHz
- Adjustable, thermostatically controlled heaters (ambient to 80°C)

SPECIFICATIONS

- Capacity: 3.25 Gallons
- Dimensions: 11.75" x 9.5" x 8"
- All Stainless Steel Enclosures
- Timer: 0 - 30 Minutes
- Lifetime Heater Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)