

DuraSonic 5.3 Gal Digital Ultrasonic Cleaner | DS20L



Stock Photo - Options May Vary

FACTORY NEW

DuraSonic DS20L 5.3 Gal Digital Ultrasonic Cleaner. Industrial Grade Model with Sweep and Degas Function and Temperature Controlled Heat for best cleaning performance. Features Operator Temperature Controls (knob), Digital Ultrasonic Timer and comes with Free Basket and Cover. Do you have questions?

To view a wide variety of Autoclaves & Equipment, click the following link to browse from automatic autoclaves and distillers to manual autoclaves and ultrasonic cleaners. Venture Medical's Complete Catalog of Autoclaves.

FEATURES

- Industrial Grade Model
- Sweep and Degas Function
- Temperature Controlled Heat for best cleaning performance
- Operator Temperature Controls (knob). Digital Ultrasonic Timer
- Comes with Free Basket and Cover

SPECIFICATIONS

- Inside Dimensions (L x W x H): 12.5" x 11" x 10"
- Outside Dimensions (L x W x H): 18"x 11" x 10.75"
- Includes Heater: Yes
- Sonic Power Watt Average: 500

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)