

DuraSonic 6.6 Gal Digital Ultrasonic Cleaner | DS25L



Stock Photo - Options May Vary

FACTORY NEW

DuraSonic DS25L 6.6 Gal Digital Ultrasonic Cleaner. The DuraSonic Digital Ultrasonic Cleaner has a 6.6-gallon (25-liter) tank, measuring 18.5" x 11.5" x 8" (L x W x D). Additionally, it includes a basket and cover, plus a 2-year warranty. Made in the USA.

Furthermore, our in-house team of certified technicians can help with any repairs or services on your current ultrasonic cleaner or autoclave. Also available are Ultrasonic Cleaner Solutions like EmPower® Dual-Enzymatic Detergent, Tuttnauer Clean & Simple Ultrasonic Enzymatic Cleaning Tablets, and SC-44 Cleaning Concentrate. TIP: Use EmPower® Foam™ Enzymatic Spray to prevent blood, tissue, mucus, and proteins from drying on instruments.

Do you have questions?

FEATURES

- Industrial Grade Model
- Sweep and Degas Function
- Digital Ultrasonic
- Comes with Free Basket and Cover

SPECIFICATIONS

- Inside Dimensions (L x W x H): 18.5" x 11.5" x 8"
- Includes Heater: Yes
- Sonic Power Watt Average: 1000

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)