

# Edan F9 Express “Loaded” Fetal/Maternal Monitor



Stock Photo - Options May Vary

## FACTORY NEW

Edan F9 Express “Loaded” Fetal/Maternal Monitor. F9 is a high-end fetal & maternal monitor; it offers the most advanced integrated monitoring for the delivery room in a big hospital, private obstetrician’s office, and antepartum clinic. F9, designed for moving situations, covers the entire continuum of antepartum, intrapartum, and postpartum applications.

The Edan F9 Express Fetal/Maternal Monitor offers a range of features for enhanced monitoring. It includes basic parameters such as FHR, TOCO, Event Mark, AFM, and internal parameters like IUP and optional DECG. Additionally, it supports 24-hour waveform playback.

You can also choose the optional built-in wireless module and software for data transmission to PC.

## FEATURES

- Basic parameters: IUP, DECG (optional)
- 4-hour waveforms playback
- Optional Built-in wireless module
- Software for data transmission to PC
- 150/152 mm wide paper

## SPECIFICATIONS

- Dimensions: 347 mm x 330 mm x 126 mm
- Weight: Approximately 6 kg
- Pulse Repetition Rate: 2KHz
- TOCO Range: 0% - 100%
- Zero Mode; Automatic/Manual

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)