

Edan SE-12 Express Stress ECG | SE-12ES



Stock Photo - Options May Vary

FACTORY NEW

Edan SE-12ES Express Stress ECG. Given the high demand within the United States, EdanUSA was formed in 2009 to offer local service and support for U.S. sales.

Additionally, all of our products come with a three-year manufacturer's warranty. We offer high-quality diagnostic equipment at 30%-40% lower prices than major U.S. competitors, improving patient care quality.

Edan SE-12ES Express Stress ECG Features: The Edan SE-12ES Express Stress ECG offers excellent performance with a built-in high-resolution thermal printer and 12 ECG traces on A4/US Letter paper. It features an 8.4" color touch-screen, alphanumeric keyboard, and one-touch operation.

FEATURES

- Excellent cost performance with Real-time waveform freezing
- Built-in high-resolution thermal printer
- Record 12 ECG traces on A4 size and US Letter-size paper
- 8.4" color touch-screen
- Alphanumeric keyboard and one-touch operation

SPECIFICATIONS

- For Complete Specifications, Please Refer to the Brochure
- Edan SE-12 Express Stress ECG Model SE-12ES Specifications (Click to View)
-

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)