

Edan SE-601C 6-Channel ECG



Stock Photo - Options May Vary

FACTORY NEW

Edan SE-601C 6-Channel ECG. Given the high demand within the United States, EdanUSA was founded in 2009 to offer local service and support for U.S. sales.

All of our products come with a three-year manufacturer's warranty. Our goal is to offer high-quality diagnostic equipment that improves the quality of patient care, with customer pricing 30% x 40% less than that of major U.S. competitors.

The Edan SE-601C 6-Channel ECG is a portable heart monitor that tracks heart rate and rhythm. It can measure heart rate variability (HRV), show R to R trends, and detect pacemakers. It has an easy-to-use alphanumeric keyboard, one-touch operation, and a rechargeable battery.

FEATURES

- Cost-effective solutions for hospital wards, routine checks, home visits, and emergency cases
- FDA Approved
- Heart rate variability (HRV) & R to R Trending analysis
- Portable, lightweight design
- Easy data input and operation

SPECIFICATIONS

- For Complete Specifications & Brochure
- Edan SE-601C 6-Channel ECG: SE-601C ([click to view](#))

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)