

Enthermics 12-Blanket Warming Cabinet | DC350



Stock Photo - Options May Vary

FACTORY NEW

The Enthermics DC350 is a 12-blanket warming cabinet designed for healthcare environments including labor and delivery rooms, ICUs, and NICUs. Its countertop design operates with no moving parts, providing quiet, non-disruptive warmth that allows patients to rest during recovery. The fan-less technology prevents lint buildup, making maintenance simple and reducing workload for healthcare staff.

The DC350 features WarmSafe™ zone heating technology that rapidly and precisely heats blankets to ideal temperatures, with heated shelves that adjust temperatures evenly throughout the cabinet. LED interior lighting with adjustable blue intensity settings creates a calming environment while allowing staff to observe inventory through the energy-efficient glass door.

The unit includes a hands-free push-button door for convenient loading and a separate warming shut-off system that prevents overheating. All chamber surface temperatures are continuously monitored to ensure efficient heat distribution and minimal energy consumption.

FEATURES

- LED interior lighting casts a comforting blue glow with two (2) different intensity settings and off-mode
- A warming shut-off system, separate from the electronic control, prevents overheating
- Glass Door: Single pane, energy-efficient, e-coated glass window in the door allows for inventory observation and reflects heat into the cavity
- Convenient: The hands-free, push-button door is designed to make loading easier when carrying an armful of blankets
- WarmRight® incorporates a multiple-zone warming technology that heats where and when it is needed
- All chamber surface temperatures are monitored, providing an efficient balance of heat, low energy consumption, and minimal heat loss
- Includes a heated center shelf to optimize heat distribution throughout the cavity further

SPECIFICATIONS

- Capacity: 10-12 blankets (3.5 ft³)
- Dimensions: 27.9 in H x 18.5 in W x 24.8 in D
- Adjustable temperature range of 37 - 82°C (98 - 180°F)
- Recommended environmental temperature range is 15°C to 32°C (60°F to 90°F)
- Relative humidity range of 10% to 95%, non-condensation
- Window: Single pane, energy-efficient, e-coated glass
- Interior/Exterior: White epoxy-coated steel exterior casing and interior insert
- Lighting: LED

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com