

Enthermics 18-Blanket Warming Cabinet | DC400



Stock Photo - Options May Vary

FACTORY NEW

The Enthermics DC400 Blanket Warming Cabinet offers a quiet, efficient solution for healthcare environments. With no moving parts, the countertop unit is perfect for labor and delivery rooms, ICUs, and NICUs. Its patient-friendly design lets you easily place it where needed, helping nurses and aides focus on care while maximizing efficiency.

The DC400 comes with several features that enhance patient comfort. Relaxing blue lights inside the unit create a calming atmosphere, contributing to a peaceful environment. WarmSafe™ zone heating technology ensures that blankets are heated quickly and precisely to the ideal temperature.

Additionally, the heated shelves maintain consistent warmth throughout the cabinet, keeping every blanket evenly warmed.

FEATURES

- LED interior lighting casts a comforting blue glow with two (2) different intensity settings and off-mode
- In the event of a power failure, the cabinet will remember its programming and begin to operate as before when power is restored
- A warming shut-off system, separate from the electronic control, prevents overheating
- The access point and the plug have been added to the back panel to allow the appliance owner to add their data logger or temperature management hardware
- Single-chamber warming cabinet White epoxy-coated steel exterior casing and interior insert

SPECIFICATIONS

- Chamber: Single-chamber warming cabinet
- Lighting: LED
- Capacity: 14-18 blankets
- Interior size: 4 ft3 (.11 m3)
- Dimensions: 21.9 in x H; 24 in W x 27 in D

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)