

# Enthermics 25-Blanket Warming Cabinet | DC750



Stock Photo - Options May Vary

## FACTORY NEW

The Enthermics DC750 25-Blanket Warming Cabinet offers a quiet, efficient solution for healthcare settings. Operating with no moving parts, this countertop unit is ideal for labor and delivery rooms, ICUs, and NICUs. Its patient-friendly design allows for easy placement in essential areas, enabling nurses and aides to focus on care while enhancing efficiency.

The DC750 is equipped with several features to improve patient comfort and care. Inside, relaxing blue lights help create a calming atmosphere, softening the room's environment. WarmSafe™ zone heating technology quickly heats blankets to precise temperatures, providing consistent warmth.

Additionally, the heated shelves adjust the temperature across the entire cabinet, ensuring each blanket warms evenly.

## FEATURES

- LED interior lighting casts a comforting blue glow with two (2) different intensity settings, and off-mode
- In the event of a power failure, the cabinet will remember its programming and begin to operate as before when power is restored. A warming shut-off system, separate from the electronic control, prevents overheating
- The access point and the plug have been added in the back panel to allow the appliance owner to add their data logger or temperature management hardware
- Single-chamber warming cabinet White epoxy-coated steel exterior casing and interior insert
- Glass Door: Single pane, energy-efficient, e-coated glass window in the door allows for inventory observation and reflects heat into the cavity

## SPECIFICATIONS

- Chamber: Single-chamber warming cabinet
- Lighting: LED
- Capacity: 20-25 blankets
- Interior size: 7.5 ft<sup>3</sup>
- Dimensions: 34.9 in H x 24 in W x 27 in D

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)