


# Mindray Passport 12 Patient Monitor with CO2 - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Mindray Passport 12 Patient Monitor with CO2 - Refurbished. The Passport 12 offers flexibility and advanced features with a modular design, making it ideal for various clinical settings.

Its plug-and-play structure allows for easy upgrades, including CO2 measurement, two additional invasive pressures, and cardiac output. Additionally, the Mindray Passport 12 features: For anesthesia monitoring, the optional multi-gas module automatically detects and measures O<sub>2</sub>, CO<sub>2</sub>, N<sub>2</sub>O, and five anesthetic agents.

Furthermore, with integrated HL7 and DIAP communication protocols, plus wired and wireless networking, it ensures smooth data sharing. The Passport 12 also connects to the BeneVision Central Station, providing centralized patient data management for better workflow and decision-making.

## FEATURES

- 12" high resolution touchscreen display
- Up to 8 waveforms standard
- Standard features include 3 or 5-lead ECG, Masimo SET SpO<sub>2</sub>, 2 IBP, NIBP, respiration, temperature, integrated recorder
- Optional Nellcor OxiMax SpO<sub>2</sub> available
- Available ST and Arrhythmia analysis with atrial fibrillation detection

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)