

# Pedigo Elite Series Warming Cabinet | P-2205



Stock Photo - Options May Vary

## FACTORY NEW

Pedigo P-2205 Elite Series Warming Cabinet. Pedigo leads in designing and constructing high-quality stainless steel and chrome products. Their products are used in emergency rooms, exam rooms, operating rooms, and sterile processing departments.

After more than a quarter-century, the Pedigo name remains synonymous with quality, durability, sterility, and exceptional customer service. Pedigo Elite Series Warming Cabinets Feature: This single-chamber warming cabinet features a white epoxy-coated steel interior insert and exterior casing. Additionally, the single-pane energy-efficient, e-coated glass window door allows easy inventory observation and hands-free door release.

Moreover, the door includes a full gasket and is hinged on the right side of the unit.

## FEATURES

- Single-chamber warming cabinet
- White epoxy-coated steel interior insert and exterior casing
- Single pane energy efficient e-coated glass window door allows inventory observation.
- Easy, hands-free door release and the door has a full gasket and is hinged on the right side of the unit
- Radiant Heat Technology utilizes multiple-zone heating that monitors surfaces for safe temperature ranges

## SPECIFICATIONS

- Temperature Range: 98°-180°
- Compartment Capacity: 1.5 cubic feet
- Compartment Dimensions: 14½"W x 15½"D x 9-5/8"H
- Overall Dimensions: 18½"W x 19-13/16"D x 16-7/8"H
- Single-chamber warming cabinet

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)