

GE/Critikon Dinamap Pro 300 Vital Signs Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

GE/Critikon Dinamap Pro 300 Vital Signs Monitor – Refurbished. The GE/Critikon Dinamap Pro 300 Vital Signs Monitor delivers dependable patient monitoring for hospitals, clinics, and outpatient facilities. Additionally, this compact monitor measures non-invasive blood pressure, oxygen saturation, and pulse rate with trusted DINAMAP technology.

Furthermore, clinicians get fast, accurate readings during routine assessments. Also, this helps medical teams maintain efficient workflows while capturing patient data throughout every shift. Moreover, the Pro 300 features an easy-to-use interface that simplifies operation in busy healthcare settings.

As a result, care teams press a single control to capture vital measurements quickly. Specifically, the monitor stores multiple readings, which helps users track patient trends during exams.

FEATURES

- Measures NIBP
- SpO2 & Heart Rate
- Pole Mount Included
- Includes: New Battery, SpO2 Sensor, NIBP Hose, 2 New Cuffs, Power Cord and Manual
- Completely Refurbished

SPECIFICATIONS

- Dimensions (without temperature): 9.7" H x 8.6" W x 5.3" D
- Dimensions (with temperature): 9.7" H x 10.0" W x 5.3" D
- Parameters: NIBP, SpO2, Heart Rate
- Included Accessories: Pole Mount, SpO2 Sensor, NIBP Hose, Two Adult Cuffs
- Battery: New

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)