

GE/Critikon Dinamap Pro 400 Vital Signs Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

GE/Critikon Dinamap Pro 400 Vital Signs Monitor – Refurbished. The refurbished GE/Critikon Dinamap Pro 400 is a reliable, cost-effective monitor that measures NIBP, SpO2, heart rate, and temperature. Additionally, this monitor is fully refurbished and biomedically inspected to ensure top performance and patient safety.

Furthermore, it works well for hospitals, clinics, and care facilities. Moreover, each Dinamap Pro 400 includes a pole mount, new battery, SpO2 sensor, and NIBP hose. Additionally, it also includes an oral temperature sensor.

As a result, staff get everything needed to start using this monitor right away. Specifically, it also comes with two adult blood pressure cuffs, a hospital-grade power cord, and a user manual.

FEATURES

- Measures NIBP
- SpO2, Temp & Heart Rate
- Pole Mount Included
- Includes New Battery, SpO2 Sensor, NIBP Hose, Oral Temp Sensor, 2 New Cuffs, Power Cord and Manual
- Completely Refurbished

SPECIFICATIONS

- Weight: 13 lbs
- Parameters: NIBP, SpO2, Heart Rate, Temperature
- Included Accessories: Pole Mount, SpO2 Sensor, NIBP Hose, Oral Temp Sensor
- Blood Pressure Cuffs: Two Adult Cuffs Included
- Battery: New

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)