

GE MAC 1200 EKG Machine with Interpretation - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Refurbished GE Marquette MAC 1200 Patient Monitoring Solution is a reliable, user-friendly system. It offers practical features ideal for smaller healthcare facilities. Designed with flexibility in mind, the MAC 1200 is a compact 12-lead electrocardiograph capable of both ECG acquisition and analysis.

The modular configuration of the MAC 1200 allows it to function as a basic ECG machine without analysis, which is ideal for facilities with minimal diagnostic needs. It can also operate as a fully featured system with interpretive capabilities. Furthermore, the MAC 1200's adaptability and proven GE Marquette technology make it a cost-effective, efficient solution for diverse clinical settings.

Do you have questions? We have answers!

FEATURES

- Full Page Printout and Preview Screen
- Waveform display
- Full-size paper reports
- Measurements, interpretation, and memory are available in any combination
- Simultaneous 12-lead ECG acquisition and analysis

SPECIFICATIONS

- Dimensions: 14.5" W x 3.7" H x 12.6" D
- Display Type: LCD 320
- x 240 pixels - Backlit contrast adjustable
- ECG Acquisition/Interpretation: 12 Lead Simultaneous
- Frequency Response: 0.04 -150 Hz

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)