

# Mortara ELI-250 EKG Machine with Interpretation - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Refurbished Mortara ELI-250 EKG Machine with Interpretation is a cost-effective, high-performance solution for your resting and diagnostic EKG needs. Furthermore, designed for efficiency and ease of use, it's the ideal EKG solution for clinics, physicians' offices, and outpatient care settings.

Additionally, its high-resolution backlit display with adjustable contrast ensures clear waveform visibility in any lighting condition. Furthermore, with intuitive one-touch operation, clinicians quickly access ECG acquisition, rhythm printing, and data transmission, streamlining workflow and saving time.

Moreover, the ELI-250 utilizes trusted Mortara ECG technology and features simple controls, providing accurate results and easy operation at an affordable price.

## FEATURES

- Full Page Printout and Preview Screen
- Waveform display
- Full-size paper reports
- Measurements, interpretation, and memory
- 12-lead EKG acquisition and analysis

## SPECIFICATIONS

- Dimensions: 15.5" W x 4" H x 17" D
- LCD Display with Backlit contrast adjustable
- ECG Acquisition/Interpretation: 12 Lead Simultaneous
- Internal Memory Storage - up to 60 ECGs
- Frequency Response: 0.05 - 300 Hz

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)