

Stryker Renaissance Advantage Series Stretchers - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker Renaissance Advantage Series Stretcher is a versatile transport and procedural stretcher designed for emergency departments, transport services, and clinical settings. It features pneumatic assist Fowler backrest, hydraulic height adjustment, Trendelenburg and reverse Trendelenburg positioning, central brake and steering, and side rails. The stretcher supports a 500 lbs weight capacity and is fully refurbished to deliver reliable performance in fast-paced clinical environments.

Key operational features include head-up and head-down positioning with pneumatic assist, adjustable height ranging from 22 to 35.5 inches, and 8-inch casters with four-wheel locking capability. The stretcher accommodates standard widths between 26 and 29 inches, depending on model. Central steering provides easy maneuverability in tight spaces.

Each unit comes with a new pad and includes a 6-month parts warranty. Multiple Stryker models are available in this series, including the 721, 1001, 1009, 1210, 1211, 1501, 1509, 1710, 1711, and 1731. IV poles are sold separately.

FEATURES

- Trendelenburg/Rev Trendelenburg
- Head-Up/Down (Pneumatic Assist)
- Height Up/Down
- Collapsible or Glide-Away Side Rails (depends on model)
- 4 Wheel Locking
- 8" Casters
- Central Steering
- Includes: New Pad

SPECIFICATIONS

- Height Range: *22" to 35.5"
- Width: *26" to 29"
- Weight Capacity: 500 lbs
- Caster Diameter: 8"
- Head-Up/Down (Pneumatic Assist)
- 6 Months Parts Warranty
- Completely Refurbished

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)