

Stryker 735 Stretcher - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker 735 Stretcher is a multi-purpose transport and care solution designed for hospitals, clinics, and surgical centers. Built with durable construction and a compact footprint, it enables easy maneuverability in tight hallways and crowded treatment areas while maintaining reliable performance in fast-paced medical environments.

The stretcher promotes both patient safety and caregiver efficiency through its intuitive design and dependable operation. Its collapsible side rails, gas spring assist backrest, and integrated features support smooth patient handling and transport.

With fifth-wheel steering, dual-end brake and steer controls, and a 500 lb weight capacity, the Stryker 735 delivers the maneuverability and stability required for emergency, transport, and general care applications.

FEATURES

- Fifth-wheel steering
- 500 lb weight capacity
- 26" wide patient surface with 4" Mattress
- Less than 21" low height
- Integrated oxygen bottle holder
- Collapsible side rails
- Four-wheel, steel-ring brake system with dual-end activators
- Dual-end brake & steer controls

SPECIFICATIONS

- 26" Litter
- Overall length: 83"
- Overall width: 30"
- Weight capacity: 500 lbs
- High Height range: 36"
- Low Height range: 21.5"
- Patient Surface: 26" x 75.5"
- Caster Diameter: 8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)