

Stryker 737 / 738 Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker 737 and 738 Stretchers are versatile, multi-purpose stretchers designed for hospitals, clinics, and urgent care centers. Built on the durable ST104 line, these stretchers feature ergonomic design and reliable performance for daily use in emergency rooms and general patient care environments.

These stretchers combine strong construction with smooth mobility to support efficient caregiver workflow while maintaining patient comfort. Key features include fifth-wheel steering for easier maneuverability in tight spaces, collapsible side rails, and integrated storage solutions.

The 737 and 738 models deliver dependable performance with a 500 lb weight capacity and multiple positioning options including height adjustment, Fowler positioning, and Trend/Reverse Trend capability. They provide a cost-effective solution without compromising on quality or versatility.

FEATURES

- Fifth-wheel steering
- 500 lb weight capacity
- Collapsible side rails
- Integrated oxygen bottle holder
- Four-wheel, steel-ring brake system with dual-end activators
- Trend/Reverse Trend
- Four integrated IV receptacles
- 8" Casters

SPECIFICATIONS

- Weight capacity: 500 lbs
- Overall length: 83"
- Overall width: 30"
- High Height range: 36"
- Low Height range: 21.5"
- Backrest: 0° - 90°
- Trend/Reverse Trend: ±18°
- Caster Diameter: 8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)