

# Stryker 747 / 748 Stretcher – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

The Stryker 747 and 748 Stretchers are multi-purpose patient transport solutions designed for durability and ease of use in fast-paced medical environments. These models from Stryker's ST104 stretcher line are ideal for emergency departments, surgical centers, and medical clinics.

Each refurbished stretcher delivers long-lasting performance and reliability without sacrificing quality or safety. The 747 features a 26" wide patient surface, while the 748 offers a 30" wide patient surface, both with 4" mattresses. Both models support a 500 lb weight capacity.

Key operational features include fifth-wheel steering for superior maneuverability, collapsible side rails, integrated oxygen bottle holder, and height adjustment ranging from 21.5" to 36". The stretchers offer Fowler and Trend/Reverse Trend positioning, along with central steering and brake controls for convenient operation.

## FEATURES

- Fifth-wheel steering
- 500 lb weight capacity
- 26" wide patient surface with 4" Mattress (747)
- 30" wide patient surface with 4" Mattress (748)
- Less than 21.5" low height
- Integrated oxygen bottle holder
- Four-wheel, steel-ring brake system with dual-end activators
- Collapsible side rails

## SPECIFICATIONS

- 26" Litter (Patient Width)
- Overall length: 83"
- Overall width: 30"
- Weight capacity: 500 lbs
- High Height range: 36"
- Low Height range: 21.5"
- Backrest: 0° - 90°
- Caster Diameter: 8"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)