

Stryker M-Series 1005 / SM104 Stretcher - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker M-Series 1005/SM104 is a bariatric stretcher designed to safely accommodate patients up to 700 pounds. It features fifth-wheel steering for enhanced maneuverability, a power-washable design, and a 26-inch wide patient surface with a new 4-inch mattress. The stretcher achieves a low height of less than 21 inches, reducing caregiver strain during patient transfers.

Key design elements include Glideaway™ side rails that create a zero-transfer gap, an integrated utility tray, and a litter-mounted oxygen bottle holder to minimize caregiver bending. The stretcher is equipped with dual side-mounted foot control hydraulics, a four-wheel steel-ring brake system with dual-end activators, and constant descent hydraulics for controlled lowering.

The M-Series offers versatile positioning capabilities with a backrest range of 0 to 90 degrees, knee gatch at 40 degrees, and trendelenburg/reverse trendelenburg positioning at ±16 degrees. Four integrated IV receptacles and an optional integrated scale system provide additional clinical functionality.

FEATURES

- Fifth-wheel steering
- 700 lb weight capacity
- 26" wide patient surface
- New 4" Mattress
- Less than 21" low height
- Litter-mounted oxygen bottle holder
- Power-washable design
- Glideaway™ side rails with integrated push/pull handles

SPECIFICATIONS

- 26" Litter
- Overall length: 84"
- Overall width, side rails up 33.5"
- Weight capacity: 700 lbs
- Low Height range: 20.75"
- Patient Surface: 26" x 75.25"
- Litter Positioning: Backrest 0° - 90°, Knee Gatch 40°, Trend./Reverse Trend.: ±16°
- Caster Diameter: 8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)