

Stryker 1105 Prime Stretcher – Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker 1105 Prime Stretcher is a versatile, multi-purpose stretcher built for demanding clinical settings. With a 700 lb weight capacity, it accommodates bariatric patients and offers two width options (26" or 30") for flexible patient care and room configurations. The 5th wheel steering system provides smooth, controlled patient transport.

Key operational features include height adjustment, Fowler positioning, and Trend/Reverse Trend capability. The stretcher incorporates a gas spring assist backrest, integrated oxygen tank holder, and glide-away side rails for ease of use. Central steering and brakes with dual-end controls enhance staff efficiency.

The stretcher combines reliability and durability with a power-washable design, four-wheel steel-ring brake system, and constant descent hydraulics. It is suitable for hospitals, surgery centers, and urgent care facilities seeking cost-effective, dependable stretcher solutions.

FEATURES

- 5th Wheel Steering
- 700 lbs weight capacity
- 26" or 30" wide patient surface with 4" Mattress
- Integrated oxygen bottle holder
- Glide-away side rails
- Four-wheel, steel-ring brake system with dual-end activators
- Trend/Reverse Trend
- Constant descent hydraulics

SPECIFICATIONS

- Overall length: 85"
- Overall width: 34" or 38"(side rails up)
- Weight capacity: 700 lbs
- High Height range: 34"
- Low Height range: 20.75"
- Patient Surface: 26" x 75.25" or 30" x 75.25"
- Caster Diameter: 8"
- Trend./Reverse Trend.: $\pm 17^\circ$

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)