

Stryker 1115 Prime Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Stryker 1115 Prime Stretcher is a multi-purpose transport solution designed for medical facilities serving bariatric patients. It features a 700 lb weight capacity, available in two width options (26" or 30"), and incorporates a Big Wheel system for enhanced maneuverability and steering compared to traditional stretchers.

The stretcher offers comprehensive height adjustment capabilities with a range from 20.75" to 34", plus Fowler and Trend/Reverse Trend positioning movements. It includes glide-away side rails, an integrated oxygen tank holder, and a gas spring-assisted backrest for improved functionality and patient comfort.

The design emphasizes durability and ease of maintenance with a power-washable construction and thermoformed ABS base. Central steering and braking controls, along with dual-end activators, provide operators with intuitive command of the stretcher during transport and positioning.

FEATURES

- Big Wheel
- 700 lbs weight capacity
- 26" or 30" wide patient surface with 4" Mattress
- Integrated oxygen bottle holder
- Glide-away side rails
- Four-wheel, steel-ring brake system with dual-end activators
- Trend/Reverse Trend
- Power-washable design

SPECIFICATIONS

- Overall length: 85"
- Overall width: 34" or 38"(side rails up)
- Weight capacity: 700 lbs
- High Height range: 34"
- Low Height range: 20.75"
- Patient Surface: 26" x 75.25" or 30" x 75.25"
- Caster Diameter: 8"
- Trend./Reverse Trend.: $\pm 17^\circ$

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)