

Stryker M-Series 1015 / SM204 Stretcher - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker M-Series 1015/SM204 is a refurbished stretcher designed to accommodate patients up to 700 pounds with enhanced safety and comfort features. Its innovative design reduces caregiver strain through user-friendly controls and mobility options while maintaining reliable patient support across diverse clinical environments.

Key features include a power-washable frame, litter-mounted oxygen bottle holder to minimize bending, and Glideaway™ side rails that create a zero-transfer gap. The stretcher offers multiple positioning options with a backrest range of 0° to 90° and integrated utility components for efficient patient care.

The M-Series incorporates Big Wheel steering that reduces startup force by 50% and steering effort by 60%, along with dual-end brake and steer controls for enhanced maneuverability. Its durable construction and versatile design make it a reliable solution for safe patient transport in modern healthcare settings.

FEATURES

- Big Wheel steering
- 700 lb weight capacity
- 26" wide patient surface
- Litter-mounted oxygen bottle holder
- Power-washable design
- Glideaway™ side rails with integrated push/pull handles
- Four integrated IV receptacles
- Dual-end brake & steer controls

SPECIFICATIONS

- 26" Litter
- Overall length: 84"
- Overall width, side rails up 33.5"
- Weight capacity: 700 lbs
- Low Height range: 20.75"
- Patient Surface: 26" x 75.25"
- Caster Diameter: 8"
- Litter Positioning: Backrest 0° - 90°, Knee Gatch 40°, Trend./Reverse Trend.: ±16°

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)