

Stryker Prime X™ X-ray Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Refurbished Stryker Prime X X-ray Stretcher is designed for seamless imaging and patient transport across clinical settings. The ClearView dual-deck design provides 360° X-ray cassette access with high-quality imaging capabilities. Three models—1105X with fifth-wheel steering, 1115X with Big Wheel system, and 1125X with Zoom Motorized Drive—accommodate diverse mobility needs.

The stretcher supports up to 700 lbs (318 kg) for bariatric care and features Glideaway Siderails for clear lateral imaging access and compact storage. The patient surface measures 76 x 191 cm with adjustable height ranging from 59.1 to 92.7 cm (23.25" to 36.5") on standard models. Backrest positioning ranges from 0° to 90°, and Trendelenburg adjustment provides ±17° capability.

User-friendly features include dual-sided hydraulic controls, integrated IV poles, oxygen bottle holders, and ergonomic push handles. Safety is ensured with a four-wheel, center-locking brake system with dual-end activators. An optional powered lift is available for the 1105X and 1115X models to facilitate effortless height adjustments.

FEATURES

- Clearview Technology Platform
- 700 lbs weight capacity
- Four-wheel, center-locking brake system with dual end activators
- Glideaway siderails
- 4 integrated IV receptacles
- Oxygen bottle holder
- Dual, side-mounted foot control hydraulics with uni-lower pedal
- +/- Trendelenburg

SPECIFICATIONS

- Patient surface: 76 x 191 cm (30" x 75.25")
- Weight capacity: 318 kg (700 lbs)
- Height range High (1105X and 1115X) 92.7 cm (36.5")
- Height range Low (1105X and 1115X) 59.1 cm (23.25")
- Backrest 0°-90° (±5°)
- Trend./Reverse Trend. (1105X and 1115X) ±17° (±3°)
- Castor diameter: 20 cm (8")
- Siderails 36 x 147 cm (14" x 57.75")

**Brochure and images are a general overview of the product. Actual product options may vary.*