

MCM 9010 Stretcher – New



Stock Photo – Options May Vary

FACTORY NEW

The MCM 9010 Medical Stretcher is a versatile patient transport and examination platform with a robust 700 lb weight capacity. It features hydraulic elevation control pedals accessible from both sides for quick height adjustment, a sliding and pneumatically assisted backrest, and support for both Trendelenburg and reverse Trendelenburg positions. Angle indicators at the head ensure accurate setup.

The stretcher accommodates examination and recovery modes with an adjustable knee gatch for patient comfort during longer stays. Fold-away push handles improve transport control, while fold-down siderails provide secure protection with quick access. A fixed IV pole with two hooks is mounted at the head end, supplemented by four IV sockets at each corner for flexible equipment placement.

Safety and functionality are enhanced by 8-inch casters with protective covers, a central-locking brake system with pedals at both ends, and bumper wheels at each corner. Additional features include drainage bag holders, a built-in oxygen tank holder, patient storage tray, and a 4-inch Comfort Series mattress pad.

FEATURES

- 700 lbs weigh capacity
- 28-inch width patient deck
- Hydraulic elevation control pedals on both sides
- Sliding backrest
- Pneumatic assisted backrest
- Trendelenburg and reverse trendelenburg
- Fold down siderails
- 4 integrated IV sockets

SPECIFICATIONS

- Weight Capacity 700 lb
- Patient surface width 28"
- Patient surface length 77"
- Height range 23.25"- 34"
- Backrest 0°- 90°
- Trendelenburg 17°
- Reverse trendelenburg -17°
- Caster diameter 8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)