

# Hausted 462 Horizon Series Stretcher - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Hausted 462 Horizon Stretcher Series is a versatile patient transport solution featuring gas spring-assisted backrest positioning, Trendelenburg and Reverse Trendelenburg capability, and hydraulic height adjustment via foot pump. Dual side foot pedal controls enable easy operation, while side rails provide patient safety during transport.

The stretcher incorporates advanced steering and braking activation from all four corners with 8-inch Tente casters for smooth maneuverability. The powder-coated steel patient surface pairs with a seamless thermoformed base cover that includes built-in storage and oxygen tank holder, maximizing clinical utility.

Additional features include a heavy-duty corner-mounted non-marring bumper system, dual accessory wells at all four corners for flexible equipment adaptation, and a 4-inch comfort mattress. The stretcher supports a maximum weight capacity of 500 lbs and accommodates both IV and accessory placement.

## FEATURES

- Side Rails
- Dual Side Foot Pedal Hydraulic Controls for Height and Trendelenburg Adjustment
- Brake and Center Track Steering Activation from All 4 Corners
- 8" Tente® Casters
- Powder Coated Steel Patient Surface
- Built-In Oxygen Tank Holder
- Heavy Duty Corner Mounted Perimeter Non-Marring Bumper System
- Dual Accessory Wells at all 4 Corners, which provide versatility in adapting accessories

## SPECIFICATIONS

- Minimum Height: 22 in
- Maximum Height: 34.75 in
- Trendelenburg: 16°
- Overall Length: 83 in
- Fowler: Pneumatic-assist (cable-free): 0-90°
- Casters: 8 in Tente
- Maximum Weight Capacity: 500 lbs
- IV and Accessory Wells: 4 IV wells / 4 accessory wells

*\*Brochure and images are a general overview of the product. Actual product options may vary.*