

Hill-Rom TranStar Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Hill-Rom TranStar Stretcher is a versatile transport stretcher designed for use in emergency departments, PACU, procedure areas, and general hospital transport. It features a pneumatic-assist Fowler backrest for smooth adjustment, hydraulic height control for easy positioning, and both Trendelenburg and reverse Trendelenburg capabilities for patient care flexibility.

The stretcher includes central brake and steer functions to enhance safety and maneuverability during transport, along with sturdy collapsible side rails for patient security. It comes with a new 4-inch pad and offers an optional knee-gatch on most models.

This refurbished unit provides dependable performance and reliability for medical facilities seeking cost-effective transport solutions. The stretcher includes a 6-month parts warranty.

FEATURES

- Completely Refurbished
- Trendelenburg/Rev Trendelenburg
- Height Up/Down
- Collapsible Side Rails
- 4 Wheel Locking
- Steering
- O2 Tank Well
- New Pad

SPECIFICATIONS

- Overall Dimensions: 26.5" or 30.25 W x 83" L
- Height Range: 21.5" - 35"
- Trend/Reverse Max Angle: 18°
- Fowler Maximum Angle: 90°
- Weight Capacity: 500 Lbs
- Collapsible Side Rails
- 4 Wheel Locking
- 6 Months Parts Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)