

Harloff Drawer Mini Line Phlebotomy Cart Pkg w/Turn latch | 3245LPB



Stock Photo - Options May Vary

FACTORY NEW

Harloff Drawer Mini Line Phlebotomy Cart Pkg w/Turn latch: 3245LPB. New Harloff Phlebotomy Packages offers a durable mobile organizational solution for blood draws, IV starts, and other specific treatment needs. Standard Features include eye-level tilt bin storage for small items, a top workspace, and a covered trash container.

In addition, Multiple locking options, cart heights, and drawer configurations are available. Alternate Uses Include Triage Supplies, Treatment, and Patient Care. Standard features include stabilizing fully curved pontoon bumpers and 5-inch premium plate casters.

Additionally, the cabinet offers 30 vertical inches of drawer space with ball-bearing full-extension slides and a turn latch.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)