

Harloff 4 Drawer Phlebotomy Cart Package w/Turn Latch, Non-Locking | 3254LPB



Stock Photo - Options May Vary

FACTORY NEW

Harloff 4 Drawer Phlebotomy Cart Pkg w/Turn Latch, Non-Locking: 3254LPB. Harloff's Phlebotomy Cart Packages offers a durable mobile organizational solution for blood draws, I.V. Starts, and other specific treatment needs.

Standard Features Include a Painted steel cabinet and drawers, eye-level tilt bin storage for small items, a top workspace, and a covered trash container. Many Drawer Configurations are Available. Alternate uses include Triage Supplies, Treatment, and Patient Care. Standard Features: Stabilizing Pontoon Bumper, 5-inch (12.7 cm) premium plate casters, Cabinet with 30 vertical inches (76.2 cm) drawer space, and Turn the latch to secure drawers while moving the cart.

FEATURES

- Stabilizing Pontoon Bumper
- 5-inch premium plate casters
- Cabinet with 30 vertical inches of drawer space
- Turn the latch to secure drawers while moving the cart

SPECIFICATIONS

- Assembled Cart Dimensions: 58.5" H x 27.5" W x 21.5" D
- Height to Work Surface - 40"
- Internal Drawer Dimensions: 15" W x 16.75" D
- Construction: Painted steel cabinet and drawers

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)