

# Harloff 5 Drawer Mini Line Phlebotomy Cart Pkg. w/Key Lock | 3245KPB



Stock Photo - Options May Vary

## FACTORY NEW

Harloff 5 Drawer Mini Line Phlebotomy Cart Pkg w/Key Lock: 3245KPB. New Harloff Phlebotomy Packages offer a durable mobile organizational solution for blood draws, I. V.

Starts, and other specific treatment needs. Standard Features include Stabilizing Pontoon Bumpers, 5" casters, key-locking system w/2 keys, and much more.

Alternate Uses Include Triage Supplies, Treatment, and Patient Care. Standard Features: Wrap-around, stabilizing bumper, 5-inch premium plate casters, Cabinet with 24 vertical inches of drawer space, Key-locking system with two keys, Pre-threaded accessory mounting holes, Ball-bearing full extension drawer slides, Removable plastic top, and Latex-free.

## FEATURES

- Stabilizing Pontoon Bumpers
- 5-inch premium plate casters
- Cabinet with 24 vertical inches of drawer space
- Key-locking system with two keys
- Pre-threaded accessory mounting holes

## SPECIFICATIONS

- Assembled Cart Dimensions: 52.3" H x 27.5" W x 21.5" D
- Construction: Painted steel cabinet and drawer
- Internal Drawer Dimensions: 15" W x 16.75" D
- Cart net weight: 108 lbs
- Approximate shipping weight: 133 pounds

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)