

# Harloff Classic 3 Drawer Isolation Cart w/Key Lock Standard Pkg | 6511



Stock Photo - Options May Vary

## FACTORY NEW

Harloff 3 Drawer Isolation Cart w/Key Lock: 6511. The Harloff 6511 Classic Three-Drawer Isolation Cart Key Lock Specialty Package delivers strong organization for isolation supplies and infection-control workflows. It supports hospitals, clinics, and treatment centers that need fast access to essential items. Furthermore, its compact design improves mobility and enhances daily efficiency.

This isolation cart includes a coat and hat rack, utility hooks, and a painted steel top rail with a push handle. The tall cabinet uses a secure key-locking system that protects critical supplies. Additionally, stabilizing pontoon bumpers increases control and shields nearby equipment.

Five-inch premium casters allow smooth transport, and the tracking and braking casters increase handling precision.

## FEATURES

- Coat and Hat Rack
- Utility Hooks
- Set-4, Stainless steel top rail/push handle
- Tall cabinet with key locking system
- The full wrap-around bumper on stabilizing the base

## SPECIFICATIONS

- Overall Dimensions: 44.5" H x 32" W x 22" D
- With coat and hat rack: 62.5" H
- Construction: All Steel
- Cart weight: 145 pounds
- Approx. shipping weight: 180 pounds

*\*Brochure and images are a general overview of the product. Actual product options may vary.*