

# Harloff Steel 6 Drawer Tall ER Crash Cart w/Emergency Package | AL809B6-EMG



Stock Photo - Options May Vary

## FACTORY NEW

Harloff 809B6-EMG Steel 6 Drawer Tall ER Crash Cart w/Emergency Package. Harloff's 20 Years of Innovation in Crash Cart Design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.

This Crash Cart comes fully equipped for demanding situations with all commonly requested accessory standards. The Harloff AL809B6-EMG features a durable powder coat paint finish over metal construction and is proudly made in the USA. Left-side pullout shelf, break-away lock bar, 100 numbered seals, 29.25" vertical drawer space, and customizable cabinet design.

Removable Plastic Top and Painted Steel Top Rail / Push Handle.

## FEATURES

- Internal Pullout shelf on the left side
- Break-away lock bar and 100 numbered break-away seals
- 29.25" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted steel top rail/push handle

## SPECIFICATIONS

- Internal Drawer Dimensions: 23" W x 17" D
- Drawer Configuration: Four-3.25", One-6.5", and One-9.75"
- Construction: Steel construction
- Overall Cart Dimensions: 67.75" H x 35" W x 24.7" D
- IV Pole extends to 85.25" H

*\*Brochure and images are a general overview of the product. Actual product options may vary.*