

Harloff Chart Holder w/Rail Clips Accessory | CHARTHLDRDM



Stock Photo - Options May Vary

FACTORY NEW

Harloff Accessory - Chart Holder w/Rail Clips: CHARTHLDRDM. The Side Chart Holder for Harloff Carts helps keep workstations organized while providing storage for charts, medical records, and documents. It conceals sensitive information in crowded areas, offering up to three inches of storage space.

This accessory mounts directly to any Harloff M-Series or A-Series medical cart and replaces accessory 680417. Call for more information. Note: Cart is Not Included, Sold Separately.

Click on the following link to view a variety of Harloff Carts. Are you looking for more cart accessories or accessories for specific needs? Then, follow the links below to view various Code/Medical/Specialty/Medication Cart Accessories.

FEATURES

- Direct Mount
- Standard sand color
- Constructed from heavy-duty steel

SPECIFICATIONS

- 8" H x 12" W x 3" D (20.3 cm x 30.5 cm x 7.6 cm)
- Weighs 3.1 lbs (1.4 kg)
- Approximate shipping weight: 7 lbs (3.2 kg)
- Approximate shipping dimensions: 17" H x 13" W x 7" D (43.2 cm x 33 cm x 17.8 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)