

Harloff OptimAL All Aluminum Punch Card Medication Cart w/Key Lock | AL460PC



Stock Photo - Options May Vary

FACTORY NEW

Harloff AL460PC OptimAL Line All Aluminum Punch Card Medication Cart with Key Lock. Nurses prefer medication carts with a two-column design for their ease of use. The narcotics and liquid storage in this med cart are separated from the rest.

This allows easy access as needed—a dual-column medication cart with 480 punch card capacity. The Harloff AL460PC Features a Single-width narcotics drawer with a BEST® brand lock on the cabinet and Illinois lock on the narcotics drawer, and Ball-bearing, full-extension, self-closing drawer slides. It also has a Plastic Replaceable Top with integrated push handles, a painted pull-out shelf, and Protective Corner Bumpers.

FEATURES

- Narcotics and liquid storage in this med cart is separated from the rest of the cart
- Dual column medication cart with 480 punch card capacity, single-width narcotics drawer
- BEST® brand lock on the cabinet and Illinois lock on the narcotics drawer
- Ball-bearing, full-extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles

SPECIFICATIONS

- Internal Drawer Dimensions: Left Side: 17" D x 20.5" W, Right Side: 17" D x 6.25" W
- Right Column Drawers: 1 ea - 4" multipurpose storage drawer, one ea - 9.75" locking narcotics drawer, and two ea - 9.75" multipurpose storage drawer
- Assembled Cart Dimensions: 45.75" H x 36" W x 23.5" D
- Writing Surface: 44" H
- Cart weight: 140 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)