

Harloff Deluxe 8 Drawer Painted Cast Cart Pkg | 6035

FACTORY NEW

Harloff Deluxe 8 Drawer Painted Cast Cart Pkg: 6035. This orthopedic cast cart system offers a convenient and cost-effective way of storing instruments and supplies for the treatment of strains, sprains, and fractures.

Designers worked with orthopedic technicians to maximize functionality.

Furthermore, The construction includes painted steel cabinets, drawers, and stainless steel accessory items.



Stock Photo - Options May Vary

FEATURES

- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with a vinyl mat
- Swivel tourniquet/pump pole

SPECIFICATIONS

- Assembled Cart Dimensions: Cart only: 43" H x 44" W x 21" D
 - Cart w/ tourniquet / pump pole: 76" H x 44" W x 21" D
 - Cart weight: 182 lbs
 - Approximate shipping weight of cart: 227 lbs
- *Structure and images are a general overview of the product. Actual product options may vary.

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)