

Harloff Deluxe 8 Drawer Painted Cast Cart w/Top Compartment | 6035-TC



Stock Photo - Options May Vary

FACTORY NEW

Harloff Deluxe 8 Drawer Painted Cast Cart w/Top Compartment: 6035-TC. This Deluxe Painted Eight Drawer Cast Cart Package with Top Compartment 6035-TC.

This orthopedic cast cart system offers a great solution to store instruments and supplies for treating strains, sprains, and fractures. Additionally, It features a painted steel cabinet, aluminum drawers, and a painted aluminum top compartment.

The Harloff 6035-TC features a top compartment for 3M® and BSN® splinting supplies with a plexiglass front. It includes a stainless steel folding side shelf, a side storage cabinet, and a scissors cable (scissors not included).

FEATURES

- Top compartment for 3M® and BSN® splinting supplies; includes plexiglass front
- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- *Scissors cable
- Stainless steel bucket holder with a vinyl mat

SPECIFICATIONS

- Cart only: 43" H x 44" W x 21" D
- Cart w/ top compartment (TC): 66" H x 44" W x 21" D
- Cart w/ tourniquet / pump pole: 76" H x 44" W x 21" D
- Cart weight: 182 lbs
- Cart weight w/ TC: 199 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)