

Harloff Deluxe 8 Drawer Stainless Steel Cast Cart Pkg | 6025



Stock Photo - Options May Vary

FACTORY NEW

Harloff Deluxe 8 Drawer Stainless Steel Cast Cart Pkg: 6025. This orthopedic cast cart system offers a convenient and cost-effective way of storing instruments and supplies for the treatment of strains, sprains, and fractures—stainless steel cabinets and drawers. The cart includes a stainless steel folding side shelf, stainless steel side storage cabinet, and scissors cable (scissors optional).

Additional 6025 Features: stainless steel bucket holder with vinyl mat, swivel tourniquet/pump pole, and 24" stainless steel ruler. The cart also includes an aluminum towel bar, internal drawer dividers, protective corner bumpers, and 5" ball-bearing casters.

One caster tracks, and one includes a brake. Full-extension ball-bearing drawer slides provide smooth operation.

FEATURES

- Stainless steel folding side shelf
- Stainless steel side storage cabinet
- Scissors cable (scissors not included)
- Stainless steel bucket holder with a vinyl mat
- Swivel tourniquet/pump pole

SPECIFICATIONS

- Assembled Cart Dimensions: Cart only: 43" H x 44" W x 21" D
 - Cart w/ tourniquet / pump pole: 76" H x 44" W x 21" D
 - Cart weight: 209 lbs
 - Approximate shipping weight: 250 lbs
- *Brochure and images are a general overview of the product. Actual product options may vary.