

Harloff Deluxe Stainless Steel 8 Drawer Cast Cart w/ Accessory Pkg & Top Storage | 6025-TC



Stock Photo - Options May Vary

FACTORY NEW

Harloff Deluxe Stainless Steel 8 Drawer Cast Cart w/ Accessory Pkg & Top Storage: 6025-TC. This orthopedic cast cart system offers a convenient and cost-effective way of storing instruments and supplies for the treatment of strains, sprains, and fractures.

Furthermore, The top compartment for 3M® and BSN® splinting supplies includes a Plexiglas front, folding shelf, and side storage cabinet. The 6025-TC features a scissors cable (scissors not included), a stainless steel bucket holder with a vinyl mat, and a swivel tourniquet/pump pole.

Additionally, it includes a 24" stainless steel ruler, an aluminum towel bar, and protective corner bumpers.

FEATURES

- Stainless steel folding shelf
- Stainless steel side storage cabinet
- Scissors cable scissors not included; see options
- Stainless steel bucket holder with a vinyl mat
- Swivel tourniquet/pump pole

SPECIFICATIONS

- Overall Dimensions: Cart only: 43" H x 44" W x 21" D
- Cart w/top compartment: (TC): 66" H x 44" W x 21" D
- Cart weight: 222 lbs
- Cart weight w/ TC: 239 lbs
- Approx *shipping weight of cart: 250 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*