

Harloff Malignant Hyperthermia Cart w/Medical Grade Refrigerator | MH5200B



Stock Photo - Options May Vary

FACTORY NEW

Harloff Malignant Hyperthermia Cart w/Medical Grade Refrigerator: MH5200B. Uses high-performance Follett® REF2 refrigerator or equivalent. Drawers allow for storing posters and MH Procedure Manual available from MHAUS.

Also, store supplies to actively cool a patient, such as hypothermia blankets, over and under the patient. A high-performance refrigerator ensures temperature controls for cold isotonic saline for IV infusions and gastric, peritoneal, or rectal irrigation. Additional MH5200B Features: The Cart can hold critical MH emergency supplies at the proper temperature in close proximity to areas where they are most needed.

This Cart has a hospital-grade plug and cord.

FEATURES

- Uses high-performance Follett® REF2 refrigerator or equivalent
- Drawers allow for storage of posters and MH Procedure Manual available from MHAUS
- Stores supplies to actively cool a patient, such as hypothermia blankets, over and under the patient
- High-performance refrigerator ensures temperature controls for cold isotonic saline for IV infusions and or gastric, peritoneal, or rectal irrigation, as appropriate
- The cart can hold critical MH emergency supplies at the proper temperature in close proximity to areas where they are most needed

SPECIFICATIONS

- Drawer Configuration: (2) - 3" and (1) - Compartment with pre-installed refrigerator
- Internal Drawer Dimensions: 23" W x 17" D
- Overall Cart Dimensions: 48.5 H x 37.5" W x 26"D
- Cart net weight: 188 pounds
- Approximate shipping weight: 227 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)