

Harloff Malignant Hyperthermia Cart w/Medical Grade Refrigerator | MH5300K



Stock Photo - Options May Vary

FACTORY NEW

Harloff Malignant Hyperthermia Cart w/Medical Grade Refrigerator: MH5300K. Uses high-performance Follett® REF2 refrigerator or equivalent.

The drawers allow for the storage of posters and the MH Procedure Manual, which are available from MHAUS. Also, store supplies to actively cool a patient, such as hypothermia blankets, over and under the patient.

Additional MH5300K Features: High-performance refrigerator ensures temperature controls for cold isotonic saline for IV infusions and gastric, peritoneal, or rectal irrigation, as appropriate. The cart can hold critical MH emergency supplies at the proper temperature in close proximity to areas where they are most needed. Breakaway locks for drawers and refrigerators with 100 breakaway seals are included.

FEATURES

- Uses high-performance Follett® REF2 refrigerator or equivalent
- Drawers allow for storage of posters and MH Procedure Manual available from MHAUS
- Stores supplies to actively cool a patient, such as hypothermia blankets, over and under the patient
- High-performance refrigerator ensures temperature controls for cold isotonic saline for IV infusions and or gastric, peritoneal, or rectal irrigation, as appropriate
- The cart can hold critical MH emergency supplies at the proper temperature in close proximity to areas where they are most needed

SPECIFICATIONS

- Construction: 20-gauge powder-coated steel with built-in refrigerator
- Overall Cart Dimensions: 51.5 H x 37.5" W x 26"D
- Cart net weight: 199 pounds
- Approximate shipping weight: 239 pounds
- Note: Refrigerator units are subject to change without notice

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)