

Harloff Malignant Hyperthermia Cart | MH5200K

FACTORY NEW

Harloff Malignant Hyperthermia Cart: MH5200K. Each Malignant Hyperthermia cart includes a built-in Follett® laboratory grade refrigerator for precise and long-lasting cold temperature control and room for required cold storage of cold saline solutions and regular insulin. Also included are drawers for storing drugs, charcoal filters, syringes, catheters, nursing supplies, and laboratory testing supplies.

The MH5200K comes with a key lock, two drawers, and a 1.8 cubic feet refrigerator. Accessories include a waste container with mounting brackets, two large chart/document holders, and a large wire basket, but an additional option is available.



Stock Photo - Options May Vary

FEATURES

- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- The highest quality in the industry, the Follett® refrigerator ensures temperature control for insulin, cold packs, and cold saline solution for IV infusions and gastric, peritoneal, or rectal irrigation, as appropriate.
- Secured drawer for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Key locks for drawers and refrigerator
- Cart ships ready to use with wheels installed

SPECIFICATIONS

- (2) - 3" (7.6 cm)
- (1) - Compartment with pre-installed refrigerator
- 23" W x 17" D (58.42 cm W x 43.18 cm D)
- Cart dims: 48.5 H x 37.5" W x 26"D (130.8 cm H x 95.3 cm W x 66 cm D)
- Cart net weight: 188 pounds (85.3 kg)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)