

Harloff Mini Wood Vinyl Al Anesthesia Cart | WV3256K-CM



Stock Photo - Options May Vary

FACTORY NEW

Harloff Mini Wood Vinyl Al Anesthesia Cart: WV3256K-CM. Medical carts with a wood furniture look are especially appealing to rehab facilities, psychiatric hospitals, and other medical facilities where the staff wants patients to feel like they're at home rather than in a sterile medical environment. This wood-finish anesthesia cart is constructed from lightweight aluminum with a coating of wood vinyl to give it a more welcoming, aesthetically pleasing look.

Features: The WV3259K Wood-Vinyl Cart has a slender build, fitting easily into narrow spaces. Its removable top adds extra writing space.

With 5" ball-bearing full swivel casters, the cart moves effortlessly through facilities.

FEATURES

- The space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Key lock with two keys
- The bumper system extends the wheelbase and stabilizes the cart

SPECIFICATIONS

- Vinyl-coated aluminum
- The top is MDF covered with a heat-pressed thermal foil laminate
- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)